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
gettyimages  
Credit: Yuichiro Chino

relationships between pupil size and certain medical conditions that affect the nervous system. Two factors, in particular, could reveal a person's mental state: the length of time it takes the pupil to dilate in response to a flash of light, and how the pupil transitions from dilation to contraction, namely whether it is progressive and smooth or halting and jagged.

For example, when a pupil does not return to its original size within six seconds after a flash of light, it's deemed not return to more than 63% of its original size, a depressive state is indicated.

Likewise, an uneven progression from contraction to dilation could indicate pain, and a rapidly alternating progression indicates the impact of medication. A post-flash dilation that exceeds the original size of the pupil can indicate excitement or nervousness.

**PUPILLOMETRY**  
Despite its potential, pupillometry is not yet part of routine medical care. Until recently, pupillometric devices cost tens of thousands of dollars and were available to few clinicians. Now, such devices are more affordable to produce and easy to transport and use.

With easier access to pupillometry, the technology could find its way into new applications. For example, AIMS is investigating the use of pupillometry as an objective measure of pain. Even the simplest discussion about pain between a patient and medical professional can be complex. Children, in particular, may not have the words to describe what they are feeling, and adults may not be able to accurately communicate the intensity of pain to their doctor, in some cases, they may

describe their pain, but not be believed. With a pupillometry device a doctor can take a rapid and reliable measure of pain, leading to either more effective treatment and reducing impact and cost. Pupillometry might be used to enhance athletic performance. Coaches can objectively measure an athlete's stress or pain levels and customize a training programme. Pupillometry devices could also be useful for tracking interest and excitement in online shopping.

**PUPILLOMETRY IS A QUICK AND NON-INVASIVE WINDOW TO THE AUTONOMIC NERVOUS SYSTEM.**

AIMS Inc. is exploring the security applications of pupillometry. Several studies have shown that the pupils of people who are lying tend to dilate more than people who are telling the truth. Used alongside other measures, they say that a pupillometry device could help detect people with criminal intentions at airports and similar public spaces that require high security.

**BRINGING AI TO THE EYE**  
Measuring how the pupil reacts is the first step, the next crucial one is analyzing the readings. To that end, AIMS has a team of five physicians and IS AI engineers and data analysts who are validating the intended applications of the company's proprietary hardware miniaturized pupillometer. To date, AIMS has received a 700-million-yen medical innovation grant from the Japanese government and applied for Pharmaceuticals and Medical Devices Agency approval for their AI-powered pupillometry device.

▲ A pupillometry device developed by AIMS Inc. in Japan.

- REFERENCES**
1. Olson, D. and Fabel, M. *Cell Case* *Ann Clin Neuro* **28**, 101-107 (2016).
  2. Hsu, L. and Hsu, M. *Science* **362**, 1491-1501 (2018).
  3. Leung, B. et al. *Present Pupillometry* **10**, 101-107 (2023).
  4. Goshima, M. et al. *Proc Natl Acad Sci* **120**, 20220101 (2023).
  5. Okubo, M. et al. *Front Res* **13**, 87-97 (2023).

www.aims-inc.jp

Advertiser retains sole responsibility for content

## FROM THE PUPIL: OF PUPILLOMETRY

KEY DEVICES THAT HARNESS AI HAVE THE POTENTIAL to make everyone

**MODERN PUPILLOMETRY**  
In 1927, Otto Lomovskoy, a professor at Columbia University in the United States, developed one of the world's first pupillometry devices, using infrared to measure the pupil's diameter. Soon after, two US-based scientists, Richard Hess and James Pitt, demonstrated that the size of the pupil changed when visual stimuli were more or less interesting, indicating a relationship between the pupil's

movement and cognitive state. Many studies since then have shown that the pupil can reflect changes in mental states, the intensity of cognition, and the nature of attention.

Over the last two decades, researchers in Japan have focused on how pupillometry can be used in diagnosis and preventative medicine," says Hiromi Kawamura, CEO of AIMS Inc. Building on this research, the healthcare startup is now investigating the

use of pupillometry as an objective measure of pain. Even the simplest discussion about pain between a patient and medical professional can be complex. Children, in particular, may not have the words to describe what they are feeling, and adults may not be able to accurately communicate the intensity of pain to their doctor, in some cases, they may

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